

1m4 Shunt band

For the secure fixation of swabs

- ▶ For securing swabs over the puncture site or wound
- ▶ Facilitates work in dialysis and clinics
- ▶ Non-slip plastic plate for optimal pressure distribution
- ▶ Adjustable elastic Velcro strap
- ▶ Skin-friendly material
- ▶ Washable at 60°C
- ▶ Custom branding available



Technical specifications

Material: Strap: 75% cotton, 20% nylon, 5% synthetic latex
Velcro: 100% nylon

Dimension: 350 x 30 mm

Weight: 14 g

HS-Code*: 63079098 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Shunt band	1 piece	25 x 1 piece	4745010564345	11000.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland