


1m4 Foot mat

To protect the padding of the examination table

- ▶ To protect against damage and contamination by shoes/feet
- ▶ Non-slip, does not slide on the bed
- ▶ Robust and durable
- ▶ Easy to clean and disinfect
- ▶ Available in four colors



Technical specifications

Material:	Polyurethane (PU), Natural rubber (NR)
Dimension:	560 x 400 x 6 mm
Weight:	640 g
HS-Code*:	40169997 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Foot mat, grey	1 piece	20 x 1 piece	4745010564413	10038.01
Foot mat, black	1 piece	20 x 1 piece	4745010564420	10038.02
Foot mat, blue	1 piece	20 x 1 piece	4745010564437	10038.03
Foot mat, burgundy	1 piece	20 x 1 piece	4745010564444	10038.04

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de