


1m4 Hand washing brush

Brush for cleaning hands and nails

- ▶ Textured handle for a secure grip
- ▶ Durable bristles
- ▶ Long-lasting
- ▶ Autoclavable, heat resistant up to 134 °C



Technical specifications

Material: Polypropylen (PP), nylon (PA)

Dimension: 105 x 36 x 35 mm

Weight: 24 g

HS-Code*: 96032980 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Hand wash. brush	1 piece	10 x 1 piece	4745010564383	10035.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de