

# 1m4 McCormick rattle

For conducting distraction hearing tests in infants

- ▶ High-frequency rattle according to McCormick
- ▶ For hearing measurement using distraction test according to Ewing
- ▶ For hearing screening in babies from 6 months as part of preventive check-ups
- ▶ Robust plastic
- ▶ Easy to clean
- ▶ Plastic free packaging
- ▶ Custom branding available



## Technical specifications

Material: Polystyrene (PS)

Dimension: 145 - Ø 60 mm

Weight: 25 g

HS-Code\*: 92089000 CN

## Product data

Product	Retail packaging	Outer packaging	EAN	REF-Nr.
McCormick rattle	1 piece	25 x 1 piece	4745010564062	10006.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: [info@medichem-online.de](mailto:info@medichem-online.de)



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland