

1m4 Instrument Brush

For cleaning instruments

- ▶ For effective cleaning of instruments
- ▶ Also removes persistent dirt
- ▶ With two nylon brush heads
Large brush head: 37 mm, bristle length: 12 mm
Small brush head: 24 mm, bristle length: 6 mm
- ▶ With ergonomic handle
- ▶ Autoclavable, heat resistant up to 134 °C



Technical specifications

Material:	Polypropylene (PP), Nylon
Dimension:	180 mm, bristle length: 6 or 13 mm
Weight:	8 g
HS-Code*:	96039099 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Double-headed Brush	1 piece	10 x 1 piece	4745010564482	12003.01

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany
Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de