


1m4 Microfilament

Diagnostic instrument for conducting sensitivity tests

- ▶ For detecting nerve damage, such as diabetic foot issues
- ▶ Protects the filament from damages
- ▶ Includes 5 replacement filaments
- ▶ Made of durable plastic
- ▶ Easy to clean
- ▶ Custom branding available



Technical specifications

Material:	Acrylonitrile-Butadiene-Styrene copolymer (ABS), Polyvinyl chloride (PVC)
Dimension:	130 x 24 mm
Weight:	15 g
HS-Code*:	90189084 CN

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Microfilament	1 piece	50 x 1 piece	4745010564123	10012.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de



TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland