

1m4 Pentherm

Diagnostic instrument for temperature sensation testing

- ▶ For early detection of distal symmetric polyneuropathy
- ▶ Based on the different thermal conductivity of the materials used
- ▶ No battery required
- ▶ Easy to clean
- ▶ With clip for fixation
- ▶ Plastic free packaging
- ▶ Custom branding available



Technical specifications

Material: Polypropylene (PP), stainless steel

Dimension: 100 x Ø 15 mm

Weight: 46 g

HS-Code*: 85131000 PK

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Pentherm	1 piece	10 x 1 piece	4745010564239	10021.00

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany

Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de

CH REP

TAS SAT AG • Chamerstr. 172 • 6300 Zug • Switzerland